

**puritymonitor**

**Version 0.1.0 (2024-09-11) unreleased**

**User manual rev. 1 (2024-09-11) unreleased**

**MIT License**

**Bastien Voirin**

**Summer 2024**

# Contents

1. Introduction and motivation .....	3
2. Conventions .....	3
2.1. Naming .....	3
2.2. Python type hints or annotations .....	3
3. Package scope .....	3
4. Package installation .....	4
4.1. Installation from PyPI or conda-forge (not implemented) .....	4
4.2. Development (editable) installation .....	4
5. Purity monitor Monte Carlo simulation .....	5
6. Purity monitor .....	6
6.1. <code>class PurityMonitor</code> .....	6
6.1.1. Constructor .....	6
6.1.2. <code>__str__()</code> .....	6
6.1.3. <code>__repr__()</code> .....	6
6.1.4. <code>draw(ax)</code> .....	6
6.2. <code>class PurityMonitorInitDecay(PurityMonitor)</code> .....	7
6.3. <code>class PurityMonitorInitDecayTimed(PurityMonitor)</code> (not implemented) .....	7
6.4. <code>class PurityMonitorFullDecay(PurityMonitor)</code> (not implemented) .....	7
6.5. <code>class PurityMonitorFullDecayTimed(PurityMonitor)</code> (not implemented) .....	7
7. Time projection chamber geometry .....	8
7.1. <code>class Geometry</code> .....	8
7.1.1. <code>__str__()</code> .....	8
7.1.2. <code>__repr__()</code> .....	8
7.1.3. <code>draw(ax)</code> .....	8
7.2. <code>class CylinderConcentricTwoPartAnode(Geometry)</code> .....	8
7.2.1. Constructor .....	8
8. TPC medium .....	9
8.1. <code>class LAr</code> .....	9
8.1.1. Constructor .....	9
9. Radioactive electron source .....	9
9.1. <code>class RadioactiveElectronSource</code> .....	9
9.1.1. Constructor .....	9
9.1.2. <code>decay()</code> .....	9
9.2. <code>class Bi207InLAr(RadioactiveElectronSource)</code> .....	9
9.2.1. Constructor .....	9
9.2.2. <code>decay()</code> .....	9
10. Command-line interface .....	10
10.1. Data file format unification .....	10
10.1.1. Data file format description .....	10
10.1.2. Data file format conversion from the oscilloscope used at B182 at CERN .....	10
10.1.3. Data file format conversion from the multichannel analyzer used at B182 at CERN .....	11
10.2. Main command .....	11
10.2.1. Macro (generic) command .....	11
10.2.2. <code>CylinderConcentricTwoPartAnode</code> geometry .....	11
11. Package structure summary .....	12
12. Example .....	12
13. Acknowledgments .....	12

# 1. Introduction and motivation

Noble liquids such as liquid argon (LAr) and liquid xenon (LXe) are chemically inert, dense, scintillating, and transparent to their own scintillation light and ionization electrons.

Successful operation of LArTPCs experiments with long drift lengths (e.g. [DUNE](#)) relies on extremely low concentrations of electronegative impurities O<sub>2</sub>, H<sub>2</sub>O, CO<sub>2</sub>, N<sub>2</sub>O...) as they hinder the free drifting of ionization electrons in the LAr volume.

## 2. Conventions

### 2.1. Naming

- `puritymonitor` (no dash, no underscore)
- PascalCase (a.k.a. UpperCamelCase) for **classes** (including `Exceptions` with the suffix “Error”) and `types`
- SCREAMING\_SNAKE\_CASE for constants
- mixedCase (a.k.a. lowerCamelCase or medial capitals) for everything else, contrary to the [PEP 8 – Style Guide for Python Code](#)

### 2.2. Python type hints or annotations

1

## 3. Package scope

The `puritymonitor` package was initially intended for R&D about a <sup>207</sup>Bi LAr PM with cylindrical TPC geometry with concentric inner disk and outer ring anodes. However, the package was designed with the intent of maximum orthogonality, including

- other TPC geometries than a cylinder with concentric inner disk anode and outer ring anode, e.g. cuboids;
- other radioactive electron sources than <sup>207</sup>Bi;
- other noble liquids than liquid argon, e.g. liquid xenon;
- other Monte Carlo simulation algorithms.

To achieve this extensibility, abstract base classes for TPC geometries, radioactive electron sources, noble liquids, or Monte Carlo simulation provide a template for derived classes to implement a particular setup.

---

<sup>1</sup>See [https://mypy.readthedocs.io/en/stable/cheat\\_sheet\\_py3.html](https://mypy.readthedocs.io/en/stable/cheat_sheet_py3.html).

## 4. Package installation

Once installed, the `puritymonitor` package can be imported in Python files using

```
import puritymonitor
```

and its command line interface (CLI) can be called from any terminal using

```
python -m puritymonitor
```

```
python -m puritymonitor -h # help
```

### 4.1. Installation from PyPI or conda-forge (not implemented)

The `puritymonitor` package can be installed from PyPI using either `pip`

```
pip install puritymonitor
```

or `conda`

```
conda install puritymonitor
```

### 4.2. Development (editable) installation

The `puritymonitor` package can be installed such as to edit its source code using `git`

```
# Clone the package repository locally without the entire git history thanks to  
# --depth 1  
git clone --depth 1 https://github.com/bastienvoirin/puritymonitor.git
```

```
# Enter the package directory  
cd puritymonitor
```

and `pip`

```
# Install the package found in the current directory in editable mode using -e  
pip install -e .
```

## **5. Purity monitor Monte Carlo simulation**

## 6. Purity monitor

### 6.1. `class PurityMonitor`

Abstract base `class`. As of version 0.1.0 of this package,

- `class PurityMonitorInitDecay(PurityMonitor)` is implemented,
- `class PurityMonitorInitDecayTimed(PurityMonitor)` is **not implemented**,
- `class PurityMonitorFullDecay(PurityMonitor)` is **not implemented**,
- `class PurityMonitorFullDecayTimed(PurityMonitor)` is **not implemented**.

However, some indications toward an implementation of the last three Monte Carlo simulations in future versions are provided in the corresponding sections.

#### 6.1.1. Constructor

```
PurityMonitor(  
    radioactiveElectronSource: RadioactiveElectronSource,  
    geometry: Geometry  
)
```

- `radioactiveElectronSource: RadioactiveElectronSource`  
Radioactive electron source (see [Section 9](#)).
- `geometry: Geometry`  
Purity monitor TPC geometry (see [Section 7](#)).

#### Example

```
purityMonitor = PurityMonitor(  
    radioactiveElectronSource = Bi207InLAR()  
    geometry = CylinderConcentricTwoPartAnode(  
        innerRadius = 10,  
        outerRadius = 20,  
        driftLength = 30  
    )  
)
```

#### 6.1.2. `__str__()`

Readable representation of the `PurityMonitor` instance.

```
str(purityMonitor)
```

#### 6.1.3. `__repr__()`

Unambiguous, explicit representation of the `PurityMonitor` instance.

```
purityMonitor.__repr__()
```

#### 6.1.4. `draw(ax)`

Draw the TPC geometry.

```
fig, ax = plt.subplots()  
purityMonitor.draw(ax)
```

## 6.2. `class PurityMonitorInitDecay(PurityMonitor)`

Only consider the electron coming from the initial de-excitation of the  $^{207}\text{Bi}$  nucleus to excited  $^{207}\text{Pb}$ . The Monte Carlo simulation thus consists in a single draw of gamma photon or internal conversion electron emitted by the initial  $^{207}\text{Bi}$  atom and ignores the subsequent de-excitations of the unstable  $^{207}\text{Pb}$  nucleus.

## 6.3. `class PurityMonitorInitDecayTimed(PurityMonitor)` (not implemented)

The arrival times of the electrons at the anode(s) are taken into account to be able to veto overlapping signals with respect to a given time delta threshold.

## 6.4. `class PurityMonitorFullDecay(PurityMonitor)` (not implemented)

The Monte Carlo simulation accounts for the full decay chain of  $^{207}\text{Bi}$  up to ground state  $^{207}\text{Pb}$ .

## 6.5. `class PurityMonitorFullDecayTimed(PurityMonitor)` (not implemented)

The arrival times of the electrons at the anode(s) are taken into account to be able to veto overlapping signals with respect to a given time delta threshold.

## 7. Time projection chamber geometry

### 7.1. `class Geometry`

Any TPC geometry must inherit from the `Geometry` abstract base class and implement the `draw(ax)` method, like the built-in `class CylinderConcentricTwoPartAnode(Geometry)` which can be used as an inspiration.

As of version 0.1.0 of the `puritymonitor` package, `CylinderConcentricTwoPartAnode` is the only built-in TPC geometry and is detailed below.

#### 7.1.1. `__str__()`

Readable representation of the `Geometry` instance.

```
str(geometry)
```

#### 7.1.2. `__repr__()`

Unambiguous, explicit representation of the `Geometry` instance.

```
geometry.__repr__()
```

#### 7.1.3. `draw(ax)`

Draw the TPC geometry.

```
fig, ax = plt.subplots()
geometry.draw(ax)
```

As `Geometry` itself is an abstract base `class`, attempting to call `draw(ax)` on an instance of `Geometry` instead of a class *derived* from `Geometry` will `raise` a `NotImplementedError`.

### 7.2. `class CylinderConcentricTwoPartAnode(Geometry)`

#### 7.2.1. Constructor

```
CylinderConcentricTwoPartAnode(
    innerRadius: float,
    outerRadius: float,
    driftLength: float
)
```

- `innerRadius: float`  
Radius of the inner disk anode (or inner radius of the outer ring anode) in `cm`.
- `outerRadius: float`  
Outer radius of the outer ring anode in `cm`.
- `driftLength: float`  
Maximum drift length, i.e. distance between the radioactive electron source and the anode plane, in `cm`.



## 8. TPC medium

### 8.1. `class` LAr

#### 8.1.1. Constructor

```
LAr()
```

## 9. Radioactive electron source

### 9.1. `class` RadioactiveElectronSource

Abstract base `class`.

#### 9.1.1. Constructor

```
RadioactiveElectronSource(  
    activity: float = float("NaN")  
)
```

- activity: `float = float("NaN")`  
Activity of the radioactive electron source in `Bq`. Only used in timed Monte Carlo simulations which are not implemented yet.

#### 9.1.2. `decay()`

### 9.2. `class` Bi207InLAr(`RadioactiveElectronSource`)

`class` derived from `RadioactiveElectronSource` defining the decay of  $^{207}\text{Bi}$  in LAr.

#### 9.2.1. Constructor

```
Bi207InLAr(  
    activity: float = float("NaN")  
)
```

- activity: `float = float("NaN")`  
Activity of the radioactive electron source in `Bq`. Only used in timed Monte Carlo simulations which are not implemented yet.

#### 9.2.2. `decay()`

## 10. Command-line interface

### 10.1. Data file format unification

A unified file format decouples the Monte Carlo simulation output and experimental data analysis input from the particular experimental setup used. Data exchange and collaboration as well as custom figure generation are all made easier by such a unified data file format, at the cost of converting experimental data coming from the acquisition apparatus first.

#### 10.1.1. Data file format description

The `puritymonitor` package expects to read and writes CSV files (i.e. comma-separated values over multiple lines) of electron energy spectra whose first column gives the lower energy bounds of the bins (i.e. rows) and header line consists in comma-separated labels for the energy column and for each spectrum column.

#### Example

```
Energy (MeV),Inner disk anode,Outer ring anode
0.0,0.0,0.0
0.01,0.0,0.0
0.02,0.0,0.0
0.03,0.0,0.0
0.04,0.0,0.0
...
0.98,80.0,80.0
0.99,90.0,90.0
1.00,100.0,100.0
1.01,90.0,90.0
1.02,80.0,80.0
...
1.96,0.0,0.0
1.97,0.0,0.0
1.98,0.0,0.0
1.99,0.0,0.0
2.0,0.0,0.0
```

#### 10.1.2. Data file format conversion from the oscilloscope used at B182 at CERN

```
puritymonitor.convert_b182_osc(*args)
```

```
python -m puritymonitor.convert_b182_osc *args
```

Short	Long	Type	Description
-d	--input-dir	path	Data files directory
-is	--inner-short	path	Short PM inner anode data filename
-il	--inner-long	path	Long PM inner anode data filename
-os	--outer-short	path	Short PM outer anode data filename
-ol	--outer-long	path	Long PM outer anode data filename
-s	--short	path	Output filename for short PM data
-l	--long	path	Output filename for long PM data

Table 1: Arguments to the `python -m puritymonitor.convert_b182_osc` command.

### 10.1.3. Data file format conversion from the multichannel analyzer used at B182 at CERN

```
puritymonitor.convert_b182_mca(*args)
```

```
python -m puritymonitor.convert_b182_mca *args
```

Short	Long	Type	Description
-d	--input-dir	path	Data files directory
-is	--inner-short	path	Short PM inner anode data filename
-il	--inner-long	path	Long PM inner anode data filename
-os	--outer-short	path	Short PM outer anode data filename
-ol	--outer-long	path	Long PM outer anode data filename
-s	--short	path	Output filename for short PM data
-l	--long	path	Output filename for long PM data

Table 2: Arguments to the `python -m puritymonitor.convert_b182_mca` command.

## 10.2. Main command

### 10.2.1. Macro (generic) command

#### 10.2.2. `CylinderConcentricTwoPartAnode` geometry

```
python -m puritymonitor.cctpa *args
```

Short	Long	Type	Description	Unit
-e	--events	int	Number of events	
-f	--field	float	Electric field	V
-a	--atten	float	Electron attenuation distance	cm
-l	--length	float	Drift length between cathode and anode planes	cm
-ir	--inner-radius	float	Inner disk anode radius	cm
-or	--outer-radius	float	Outer ring anode radius	cm
-r	--relative-scale	list[float]	Relative scaling between the inner disk anode spectrum and outer ring anode spectrum.	
-g	--geom	Flag	TPC geometry visualization	
-d	--data	Flag	Experimental data	
-s	--simu	Flag	Monte Carlo simulation	

Table 3: Arguments to the `python -m puritymonitor.cctpa` command.

## 11. Package structure summary

```
class PurityMonitor
    class PurityMonitorInitDecay(PurityMonitor)
    class PurityMonitorInitDecayTimed(PurityMonitor)
    class PurityMonitorFullDecay(PurityMonitor)
    class PurityMonitorFullDecayTimed(PurityMonitor)
```

```
class RadioactiveElectronSource
    class Bi207InLAr(RadioactiveElectronSource)
```

```
class Geometry
    class CylinderConcentricTwoPartAnode(Geometry)
```

```
class LAr
```

```
class EnergySpectra
```

## 12. Example

## 13. Acknowledgments

This Python package has been developed as part of the **CERN Summer Student Programme 2024** under the supervision of **Francesco Pietropaolo** ([CERN](#)) from his Monte Carlo simulation script, data analysis scripts and experimental data from **Robert Gan** ([Boston University](#)), and experimental data from **Gajendra Gurung** ([University of Texas at Arlington](#)) with the intent of the  $^{207}\text{Bi}$ -based LAr purity monitor R&D at CERN and INFN Padova.